
## Pseudalert\* Quick Reference Guide

## Quick tips for reading Pseudalert results

After 24 hours of incubation at 38°C  $\pm 0.5$ °C (or 26 hours for carbonated samples):

- Read results using a 6-watt, 365-nm ultraviolet (UV) lamp<sup>1</sup>; ensure the lamp has the proper wavelength and intensity.
- Use a viewing cabinet<sup>2</sup> or read results in a dark room to avoid interference from ambient light.
- Interpret any blue fluorescence as **POSITIVE**, even if the fluorescent signal is weak.
- If you are unsure about a well or vessel with weak fluorescence at 24 hours, incubate for another 1–4 hours.

Fluorescence Read Tips (as seen under UV light)



For more information, visit us at idexx.com/pseudalert. In North America, call 1-800-321-0207. In the United Kingdom, call +44 1638 676 800.

