HPC for Quanti-Tray*

IDEXX

06-27391-02

idexx.com/water

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

IDEXX

Australia: 1800 443 399

Japan: +81-422 71 5921

China: +86-21 61279528

UK: +44 (0) 1638 676800

Europe: 00800 4339 9111

North/South America: 1 207 556 4496/1 800 321 0207

For Technical Support, please call:

HPC for Quanti-Tray* Test Kit

Introduction

The HPC for Quanti-Tray* Test is used for the quantification of heterotrophic organisms in 100 mL of water. The HPC for Quanti-Tray Test is based on IDEXX-patented Multiple Enzyme Technology*, which detects viable heterotrophic organisms by testing for the presence of key enzymes known to be present in these organisms. The test uses multiple enzyme substrates that produce blue fluorescence when metabolized by waterborne bacteria. The sample/reagent mixture is added to a Quanti-Tray*, incubated, and then examined for fluorescing wells. The number of fluorescing wells corresponds to a most probable number (MPN) of total heterotrophic organisms in the original sample. The HPC for Quanti-Tray Test detects organisms at 1 cfu/100 mL after 44–72 hours of incubation.

Storage

Store at 2–25°C away from light.

Quanti-Tray Enumeration Procedure

1. Add contents of one pack to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until powder is dissolved.
3. Pour sample/reagent mixture into a Quanti-Tray or Quanti-Tray*/2000 and seal using a Quanti-Tray* Sealer.
4. Place the sealed tray in a 36° ± 2°C incubator for 44 hours. Results are valid for up to 72 hours.
5. Interpret results according to the Result Interpretation table below. Refer to the MPN table provided with the trays to obtain a most probable number.

Result Interpretation

Appearance	Result
Blue fluorescence [†]	Positive for heterotrophic organisms
No blue fluorescence	Negative for heterotrophic organisms

[†]Presence of blue fluorescent color greater than the amount present in a negative control sample.

- Look for fluorescence with a 6-watt, 365-nm UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- HPC for Quanti-Tray results are definitive at 44–72 hours.

Procedural Notes

- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- HPC for Quanti-Tray sample/reagent mixture will appear pale pink in color after the reagent is added to the sample. This color helps to differentiate the HPC for Quanti-Tray samples from other IDEXX water test samples in the laboratory.
- If the sample has been diluted, multiply the MPN value by the dilution factor to obtain the correct quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- The HPC for Quanti-Tray Test is a primary water test. HPC for Quanti-Tray performance characteristics do not apply to samples altered by pre-enrichment or concentration.
- Chlorinated samples should be treated with sodium thiosulfate prior to testing. If the chlorine is not neutralized, consider the sample invalid and discontinue testing.
- For comparison, an incubated sterile water blank containing HPC for Quanti-Tray reagent (negative control) can be used when interpreting results.
- Aseptic technique should always be followed when using the HPC for Quanti-Tray Test. Dispose of in accordance with Good Laboratory Practices.

Quality Control Procedure

The following quality control procedures is recommended for each lot of the HPC for Quanti-Tray Test:

1. Positive Control: IDEXX-QC HPC/TVC¹: *Enterococcus faecalis*.
 - A. Review the range for the control - if greater than 200, a dilution may be necessary for the Quanti-Tray prior to proceeding to step 2.A.
2. Negative/blank control:
 - A. Add the contents of one snap pack to 100 mL of sterile, deionized water.
 - B. Cap the vessel and shake until powder is dissolved.
 - C. Pour sterile water/reagent mixture into a Quanti-Tray or Quanti-Tray/2000, and seal using a Quanti-Tray Sealer.
 - D. Incubate the sealed tray at 36° ± 2°C for 44–72 hours.
 - E. Note appearance of wells under UV light. Negative Control/Blank wells should not fluoresce up to 72 hours of incubation.

1. IDEXX-QC HPC/TVC, IDEXX Catalog # UN3373-WQC-HPC

*Quanti-Tray and Multiple Enzyme Technology are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2013 IDEXX Laboratories, Inc. All rights reserved.

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© 2013 IDEXX Laboratories, Inc. Todos os direitos reservados.