
Enterolert^{*}-E

06-04626-09

IDEXX

For Technical Support, please call:

North/South America: 1 207 556 4496/1 800 321 0207

Europe: 00800 4339 9111

UK: +44 (0) 1638 676800

China: +86 21 61279528

Japan: +03 5301 6800

Australia: 1300 443 399

IDEXX

IDEXX Laboratories, Inc., One IDEXX Drive, Westbrook, Maine 04092 USA

idexx.com/water

Enterolert[®]-E Test Kit

Introduction

Enterolert-E detects enterococci, such as *E. faecium* and *E. faecalis*, in fresh and marine water. It is based on IDEXX's proprietary Defined Substrate Technology[®] (DST[®]). When enterococci utilize their β -glucosidase enzyme to metabolize Enterolert-E's nutrient-indicator, 4-methyl-umbelliferyl β -D-glucoside, the sample fluoresces. Enterolert-E detects enterococci at 1 cfu per 100 mL sample within 24 hours.

Storage

Store at 2–30°C away from light.

Presence/Absence (P/A) Procedure

1. Add contents of one pack to a 100 mL sample in a sterile, transparent, non fluorescing vessel.
2. Cap vessel and shake.
3. Incubate at 41 ± 0.5°C for 24 hours.
4. Read results according to Result Interpretation table below.

Quanti-Tray[®] Enumeration Procedure

1. Add contents of one pack to a 100 mL water sample in a sterile vessel.
2. Cap vessel and shake until dissolved.
3. Pour sample/reagent mixture into a Quanti-Tray[®] or Quanti-Tray[®]/2000 and seal in an IDEXX Quanti-Tray[®] Sealer.
4. Place the sealed tray in a 41 ± 0.5°C incubator for 24 hours.
5. Read results according to the Result Interpretation table below. Count the number of positive wells and refer to the MPN table provided with the trays to obtain a Most Probable Number.

Result Interpretation

Appearance	Result
Lack of fluorescence	Negative for enterococci
Blue fluorescence	Positive for enterococci

- Look for fluorescence with a 6-watt, 365 nm, UV light within 5 inches of the sample in a dark environment. Face light away from your eyes and towards the sample.
- Enterolert-E results are definitive at 24–28 hours. In addition, positives for enterococci observed before 24 hours and negatives observed after 28 hours are also valid.

Procedural Notes

- Marine water samples must be diluted at least tenfold. When following AFNOR validated method for bathing water, samples must also be diluted at least tenfold. Multiply the MPN value by the dilution factor to obtain the proper quantitative result.
- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a water blank can be used when interpreting results.
- Enterolert-E can be run in any multiple tube format. *Standard Methods for the Examination of Water and Wastewater*¹ MPN tables should be used to find Most Probable Numbers.
- Enterolert-E is a primary water test. Enterolert-E performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- Aseptic technique should always be followed when using Enterolert-E. Dispose of in accordance with Good Laboratory Practices.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.

Quality Control Procedures

1. One of the following quality control procedures is recommended for each lot of Enterolert-E:
 - A. IDEXX-QC Enterococci²: *Enterococcus faecalis*, *Escherichia coli*, *Streptococcus bovis*.
 - B.
 - i. For each of the American Type Culture Collection (ATCC)³ bacterial strain (*Enterococcus faecium* ATCC 35667, *Serratia marcescens* ATCC 43862, and *Aerococcus viridans* ATCC 10400), streak the culture onto labeled TSA or Blood Agar plates and incubate at 35 ± 2°C for 18–24 hours.
 - ii. For each bacterial strain, touch a sterile 1 μ L inoculating loop to a colony and use it to inoculate a labeled test tube containing 5 mL of sterile deionized water. Close cap and shake thoroughly.
 - iii. For each bacterial strain, take a 1 μ L loop from the test tube and use it to inoculate a labeled vessel containing 100 mL of sterile deionized water. These are your controls.
2. Follow the P/A Procedure or Quanti-Tray Enumeration Procedure above.
3. Results should match the Result Interpretation table above.

NOTE: IDEXX internal quality control testing is performed in accordance with ISO 11133:2014. Quality Control Certificates are available at idexx.com/water.

1. Eaton, AD, Clesceri, LS, Greenberg, AE, Rice, EN. *Standard Methods for the Examination of Water and Wastewater*. American Public Health Association, 2005. Washington, DC.

2. IDEXX-QC Enterococci, IDEXX Catalog #UN3373-WQC-ENT

3. American Type Culture Collection 1-800-638-6597 www.atcc.org

*Enterolert, Defined Substrate Technology, DST and Quanti-Tray are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

Patent information: idexx.com/patents.

© 2017 IDEXX Laboratories, Inc. All rights reserved.

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1. Eaton, AD, Clesceri, LS, Greenberg, AE, Rice, EN. *Standard Methods for the Examination of Water and Wastewater (Standardverfahren für die Wasser- und Abwasseruntersuchung)*. American Public Health Association, 2005. Washington, DC, USA.

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Patentinformation: idexx.com/patents.

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IDX 33/04 - 02/15
WATER ANALYSIS METHODS
<http://nf-validation.afnor.org>

The method Enterolert*-E/Quanti-Tray* or Quanti-Tray*/2000 for water analysis is granted NF Validation by AFNOR Certification as an alternative method to the standard ISO 7899-1 for enumeration of Intestinal enterococci in bathing water, under the Certificate number: IDX 33/04-02/15.

For more information about end of validity, please refer to the certificate NF Validation available on website mentioned above.

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IDEXX Water Quality Control Laboratory is accredited to ISO/IEC 17025:2005

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© 2017 IDEXX Laboratories, Inc. Alle Rechte vorbehalten.