

SNAP™ FeLV Antigen Test

Test accuracy

Sensitivity and specificity of the SNAP™ FeLV Antigen Test

Analyte	Reference method	SNAP test result		Total	Sensitivity (95% CL)
		+	–		Specificity (95% CL)
FeLV ¹	+	158	0	407	100% (97.7%–100%)
	–	2	247		99.2% (97.1-99.9%)

Source: Singh S, Davenport KA, Schooley E, Ruggiero A, Nassar S, Buch J, Chandrashekar R. Diagnostic Accuracy of a Point-of-Care Immunoassay for Feline Immunodeficiency Virus Antibodies, Feline Leukemia Virus Antigen, and *Dirofilaria immitis* Antigen. *Viruses*. 2023 Oct 19;15(10):2117. doi:10.3390/v15102117. PMID:37896893; PMCID: PMC106112086.

Reference method

1. PetChek* Feline Leukemia Virus Antigen Test

For interpretation and diagnostic algorithms, visit idexx.co.uk/FelVAntigenTestResources.