

# The accurate Catalyst SDMA Test you know and trust. Now upgraded.

**Easier to use. Faster results. And simple freezer storage with other routine slides.**

Kidneys are earlier indicators of many aspects of kidney health. By including the new Catalyst\* SDMA Test with every chemistry panel, it's now easier than ever to hear everything kidneys have to say.

Find out more about the Catalyst SDMA Test improvements, watch a short video on how easy it is to run, read FAQs and get staff training at [idexx.co.uk/newSDMA](https://idexx.co.uk/newSDMA).



## What's changed

- + Convenient storage: Now in the freezer, making it easier to grab and go when getting other Catalyst\* slides.
- + Simpler workflow and faster results: No more reagent cups. Now that the reagent has been incorporated into the slide, you can simply add the slide to your chemistry run for faster results.
- + New product numbers: Just sign in to IDEXX Online Orders to see all the ways to order SDMA with your other chemistries.
- + Ordering kits: The Catalyst\* SDMA and Total T<sub>4</sub> Kit will be discontinued, but the tests can still be purchased as a bundle.

## What's staying the same

- + The industry-leading, reliable accuracy you've come to expect so you can gain deeper insights into kidney function.
- + Ordering convenience: SDMA will still be included in your most commonly run profiles—get everything you need in one click.
- + Flexibility to add to any panel in any order.
- + Same quantity, less bulk: You'll still get 12 tests per box, but smaller packaging means easier storage and less waste.

